

# Thread Lifting



# Tentech Item

10 cog x 10 mono  
Thread Thread

<p><b>10 THERMA</b></p> <p>Monopolar RF</p>	<p><b>10 THERA</b></p> <p>2 Line HIFU</p>	<p><b>10 PL</b></p> <p>168 Cell Fractional IPL</p>	<p><b>10 CELL</b></p> <p>3cm x 3cm CO2 Fractional Laser</p>	<p><b>ONE THERA</b></p> <p>Dual Mode Ultrasound System</p>	<p><b>TEN HI</b></p> <p>Fractional &amp; Thermal RF</p>	<p><b>10 SONO</b></p> <p>DSS &amp; NBS Dual Mode System</p>	<p><b>LIPO THERA</b></p> <p>Dual Mode Ultrasound System</p>	<p><b>10 LUX</b></p> <p>High Intensity &amp; Multiple Light Therapy Device</p>	<p><b>10 LUX HAIR</b></p> <p>PhotoBioModulation Device</p>
				<p><b>Rhino fill</b> Premium</p> <p><b>Regevan fill</b> volume MAX</p> <p><b>Regevan fill</b> volume PRO</p> <p><b>Regevan fill</b> volume UP</p> <p><b>Regevan fill</b> Fine</p> <p>Monophasic Filler [for Export]</p>		<p><b>Regevan<sup>Ph</sup></b></p> <p>Polynucleotide Sodium</p>			
<p><b>Regevan<sup>LED</sup> Mask</b></p> <p>Dr. Oracle LED Mask 4,474 LED Lights</p>	<p><b>10 mono</b> Thread</p> <p>Mono Thread [Registered for NMPA]</p>	<p><b>10 cog</b> Thread</p> <p>Asymmetrical Cog Thread</p>					<p><b>Regevan<sup>C</sup></b></p> <p>Growth Factor Ointment</p>	<p><b>Regevan<sup>S</sup></b></p> <p>Growth Factor Spray</p>	<p><b>Regevan<sup>ALO CARE</sup></b></p> <p>Growth Factor Ampoule</p>

- **What is Thread lifting?**  
a procedure for lifting the skin using thread
- **Principle of Thread Lifting**  
The Thread itself has a cog, which physically pulls tissue directly, and as the Thread melts, collagen is produced around it, resulting in elasticity



## 10 mono Thread


Ingredient	Needle Length	Thread Length	USP	Package
Type	120mm	190mm	1-0	10ea
Component	80mm	150mm	2-0	10ea

## 10 cog Thread

Needle Size (Gauge)	Needle Length	Thread Length	USP	Package
29g	25mm	40mm	6-0	20ea
29g	38mm	50mm	6-0	20ea
26g	60mm	90mm	5-0	20ea


## 1) Easy procedure

Due to the manufacture of ultrasonic molding method, the thread is strong, so even if it is asymmetric, it maintains the pulling force of the tissue. Easy procedure to insert the needle in the direction of lifting from the chair position




## 2) Manufactured by ultrasonic processing

Overcoming the disadvantages of heat treatment with 3D compression using ultrasonic molds. No thermal toxicity by ultrasonic method other than normal heat treatment



**VS**



# 10 cog x 10 mono

Thread Thread

All information in this document is subject to copyrights and other intellectual property rights of Tentech Inc. and other licensees. No modification, reproduction or copy of any or any part of this material may be made without prior written permission from Tentech Inc. and licensees. Additional information is available upon request. The following are trademarks or registered trademarks of Tentech, Inc. in Korea and other countries.

